


# SaliVIR COVID19 Rapid 1 STEP QRT-PCR kit



## Chronomics Protocol Guide v1.0.2

### SaliVIR COVID19 Rapid 1 STEP QRT-PCR Protocol

#### Intended Use

The SaliVIR COVID19 Rapid 1 STEP QRT-PCR kit is intended to be used for the detection of SARS-CoV-2 genomic RNA extracted from biological samples derived from lower or upper respiratory tract specimens. The inhibitor tolerant qRT-PCR mix provided in the kit is a one-step solution designed for amplification of the SARS-CoV-2 RNA targets and does not contain an internal reference dye. Included in the kit are the N1 and N2 primer / probe assays targeting the nucleocapsid gene of SARS-CoV-2; and the RP primer / probe assay to be used as internal RNA extraction control. Two plasmid controls are also included in the kit to confirm functionality of the assays and the qRT-PCR reaction: the RPP30 Negative Control and the 2019-nCoV nucleocapsid gene Positive Control.

#### Kit Contents

Reagent	Volume (1,000 rxns)	Volume (10,000 rxns)	Storage
SaliVIR qRT-PCR Master Mix	1.5 ml	2 x 30 ml	- 20°C
N1 primer / probe mix *	2 x 0.9 ml	1 x 18 ml	- 20°C
N2 primer / probe mix *	2 x 0.9 ml	1 x 18 ml	- 20°C
RP primer / probe mix *	2 x 0.9 ml	1 x 18 ml	- 20°C
RPP30 Negative Control	1 x 0.5 ml	2 x 1 ml	- 20°C
2019-nCoV nucleocapsid gene Positive Control	1 x 0.5 ml	2 x 1 ml	- 20°C

\* All primer / probe mix assay tubes contain a 6.7 µM concentration of each primer (i.e. forward and reverse) and a 1.7 µM concentration of probe. Once added to the reaction mix, primers will be present at a working concentration of 500 nM and the probe at a working concentration of 125 nM.

#### Required equipment

- Class II Biological safety cabinet
- Single and/or multichannel pipettes (10, 100, 200, 1000 µl)
- PCR-clean filtered tips
- 1.5 / 2 ml cold block (or access to ice)
- 96 well cold block (or access to ice)
- qPCR Instrument (4 colour)
- 96 well plate and optical seal compatible with qPCR instrument

#### Additional User Supplied Consumables

- Molecular biology grade water

#### Storage and Handling

Upon receipt, store all reagents at -20°C.

Thaw the SaliVIR qRT-PCR Master Mix on ice and keep on ice at all times. Thaw the RP, N1 and N2 assay primer / probe mixes at room temperature in the dark and then keep on ice throughout the qRT-PCR setup preparation. After thawing, ensure that all reagents are mixed briefly by vortexing and then spun down. Avoid repeated freeze / thawing whenever possible.

## Handling of Plasmid Controls

The RPP30 Negative Control and 2019-nCoV nucleocapsid gene Positive control consist of plasmids containing the RPP30 gene and the nucleocapsid gene from SARS-CoV-2 respectively and are provided at a concentration of 200 copies/ $\mu$ L. When thawed for the first time, it is recommended to aliquot out the entire amount of both controls in single use aliquots in order to minimize freeze / thaw cycles. For each qRT-PCR run, 15  $\mu$ L of each control are used (i.e. 5  $\mu$ L per reaction = 1000 copies per reaction).

## qRT-PCR Setup Procedure

Table 1 outlines the required volumes of sample / control and reagents needed to set up the qRT-PCR reaction for one assay. Prepare three reaction mixes separately; one reaction mix for assay N1, one for assay N2 and one for assay RP. Each sample / control requires one replicate for each assay. Each qRT-PCR run should include the No Template Control (i.e. molecular biology grade water), the RPP30 Negative Control and the 2019-nCoV nucleocapsid gene Positive Control. All procedures should be carried out in a sterile environment, ideally a Class II biosafety cabinet. Thaw reagents as described above and setup the reaction on ice (or on cold block).

**Table 1. qRT-PCR reaction mix setup volumes for a single reaction. When preparing a mix for multiple reactions, include a 5% overage for each reagent.**

Reagent	Volume
SaliVIR qRT-PCR Master Mix	5 $\mu$ L
Primer and Probe Mix	1.5 $\mu$ L
Template	5 $\mu$ L
Water	8.5 $\mu$ L
<b>Total</b>	<b>20 <math>\mu</math>L</b>

**NOTE:** the volume of water and template can be adjusted to include more template in the reaction mix as required.

## qRT-PCR Setup and Cycling Conditions

Set up the qPCR instrument using manufacturers guidelines. Where possible, choose the Quantitation by Comparative Ct ( $\Delta \Delta$ Ct) method with TaqMan or "Other" reagents (do not add a melt curve option). For qPCR instruments with "FAST" blocks, select the Standard ramp speed. Select the FAM filter for the N1, N2 and RP targets. If possible, select NFQ-MGB as quencher; or alternatively leave this field empty. Do not select a dye as quencher (such as TAMRA). Select Program cycling conditions as shown in Table 2 below. Set reaction volume to 20  $\mu$ L.

**Table 2. qRT-PCR Program cycling conditions**

Step	Cycles	Temperature	Time
1	1	50 °C	10 min
2	1	95 °C	2 min
3	45	95 °C	5 s
4		55 °C	30 s

**NOTE:** fluorescence acquisition is performed at step 4.

**qRT-PCR Data Analysis and Interpretation**

Please note data analysis may vary between qPCR machines and thresholds must be determined empirically by the end user or laboratory. We recommend setting the Baseline start cycle at 5 and the end cycle at 15; and the threshold at 200 RFU or 0.02 ΔCt.

As per CDC guidelines, Ct values that fall below the 40 cycles threshold are considered positive signals. Refer to Table 3 below for interpretation of results from control and patient samples.

**Table 3. Interpretation of results from control and patient samples**

Sample	RP Result	N1 Result	N2 Result
RPP30 Negative Control	+	-	-
2019-nCoV nucleocapsid gene Positive Control	-	+	+
Positive patient sample	+	+	+
	-	+	+
Negative patient sample	+	-	-
Inconclusive patient sample	+	+	-
	+	-	+
	-	+	-
	-	-	+
Failed patient sample	-	-	-

**Kit Specification and Performance**


<b>Application</b>	Qualitative PCR test for detection of SARS-CoV-2 nucleocapsid gene
<b>Type of Detection</b>	Ribonucleic acid (RNA) of SARS-CoV-2
<b>Sample Type</b>	Lower respiratory tract specimens (e.g bronchoalveolar lavage, sputum, tracheal aspirate) and upper respiratory tract specimens (e.g saliva, nasopharyngeal fluids, nasal swab)
<b>qRT-PCR Limit of Detection</b>	1x10 <sup>0</sup>

### Limit of Detection (LoD) Testing


To understand the limit of detection (LoD) for the SaliVIR COVID19 Rapid 1 STEP QRT-PCR kit we undertook the following experiment to calculate experimentally the minimum copies of viral template RNA which can be detected.

Synthetic SARS-CoV-2 RNA: ORF, E, N (ATCC® VR-3276SD™) was used in a serial dilution experiment to establish a limit of detection of SaliVIR qRT-PCR Master Mix using N1 and N2 primer and probe mixes on a BioRad CFX96.

**N1 Assay**



**N2 Assay**



Sample	Copies	Type	N1 Cq	N1 Cq StDev	N2 Cq	N2 Cq StDev
1x10 <sup>0</sup>	1	Synthetic RNA	37.01	1.16	37.45	0.93
1x10 <sup>0.7</sup>	5	Synthetic RNA	35.61	0.31	36.77	1.22
1x10 <sup>1</sup>	10	Synthetic RNA	34.92	0.76	35.18	0.22
1x10 <sup>2</sup>	100	Synthetic RNA	31.00	0.31	31.46	0.32
1x10 <sup>3</sup>	1000	Synthetic RNA	27.58	0.10	28.05	0.12

